

CHAPTER I

INTRODUCTION

1.1 Background

The stirred tank reactor is the most common reactor in the chemical industry. In large-scale industries, stirred tank flow reactors are more often applied because of their adjustable operating capacity. The performance of the stirred flow reactor needs to be studied to determine the characteristics of the fluid flow, the reactions that occur in the optimization of the operation of the reactor.

The operation of the stirred tank flow reactor includes three stages, namely filling the reactor with high overflow, continuous conditions, and continuous steady state. Evaluation of operating variables is very easy to do at steady state conditions.

Mathematical modeling is needed to facilitate the analysis of problems that arise in the operation of the stirred tank flow reactor. The proposed mathematical model is tested for accuracy by comparing it with experimental data. The proposed mathematical model is solved by means of analysis if the equation is easy to solve. However, for complex reactions, a complex mathematical model will also be obtained. Numerical method is highly recommended to obtain the value of k , mass transfer constant, and reaction order which are adjustable parameters.

1.2 Problem Formulation

The Ideal Continuous Flow Reactor / Stirred Tank Flow Reactor is a place where a chemical reaction of formation or decomposition occurs, where the mass flow enters or exits continuously (continuously) (Rosadi, 2000). In this practicum, the ideal continuous flow reactor will be studied to determine the value of the reaction order, reaction constants and the effect of reactant concentration and compare the calculation of the mathematical model of ethyl acetate saponification with NaOH in an Ideal Continuous Flow Reactor.


1.3 Purpose

1. Determine the value of the reaction order for ethyl acetate saponification with NaOH.
2. Calculate the value of the reaction constant (k) for saponification of ethyl acetate with NaOH.

3. Knowing the effect of variables on the reaction constant (k) of saponification of ethyl acetate with NaOH.
4. Comparing the experimental results with the calculation of the mathematical model of the Runge Kutta method for the saponification reaction in an ideal continuous flow reactor.

1.4 Benefits

1. Students can determine the order of the reaction for saponification of Ethyl acetate with NaOH.
2. Students can calculate the value of the reaction constant (k) for saponification of Ethyl acetate with NaOH.
3. Students are able to determine the effect of variables on the reaction constant (k) of saponification of ethyl acetate with NaOH.
4. Students are able to compare the experimental results with the calculation of the mathematical model of the Runge Kutta method for the saponification reaction in an ideal continuous flow reactor.



Process

Laboratory

CHAPTER II

LITERATURE REVIEW

2.1 Batch Reactor

Material balance in reactor simultaneously



Figure 2.1 Part of the mass balance of a system

Input = 0

Output = 0

Reactants that react = $(-r_A)$

Input = output + reactants react + accumulation

$$0 = 0 + v(-r_A) + \frac{dN_A}{dt} \quad (2.1)$$

$$0 = V_i(-r_A) + \frac{d[N_{A0}(1-X_A)]}{dt} \quad (2.2)$$

$$0 = V_i(-r_A) - \frac{N_{A0} dX_A}{dt} \quad (2.3)$$

$$dt = \frac{N_{A0}}{V_i(-r_A)} dX_A \quad (2.4)$$

$$t = N_A \int_0^{X_A} \frac{dX_A}{V_i(-r_A)} \quad (2.5)$$

At constant volume

$$\begin{aligned} C_A &= C_{A0}(1-X_A) \\ dC_A &= -C_{A0} \cdot dX_A \end{aligned} \quad (2.6)$$

Substituting Eq. (6) into Eq. (5) gives:

$$t = C_{A0} \int_0^X \frac{dN_A}{-r_A} = - \int_{C_{A0}}^{C_A} \frac{dC_A}{-r_A} \quad (2.7)$$

2.2 Continuous Flow Ideal Reactor/Stirred Flow Tank Reactor (CSTR)

The stages that occur in the CSTR reactor are divided into 3 process stages, namely:

a. The first stage

The first stage starts at $t=0$ until overflow

From the law of conservation of mass

Accumulation = Input – Output

$$\rho \frac{dV}{dt} = \rho \cdot F_0 - 0 \quad (2.8)$$

$$dV = F_0 \cdot dt, \text{ at } t = 0 \rightarrow V = 0$$

Since the density of the flow rate is considered constant, its volume is only a function of time.

$$V = F_0 \cdot t \quad (2.9)$$

Meanwhile, from the component balance sheet

Accumulation = input – output – rate of consumption due to reaction

$$\frac{d}{dt}(V \cdot C) = F_0 \cdot C_0 - 0 - V(-r_A) \quad (2.10)$$

In this case:

V = Volume material in reactor (l)

C = Molar concentration of reactants in the reactor (mol/l)

F_0 = Incoming reactant flow rate (l/mol)

C_0 = Molar concentration of reactants in feed (mol/l)

t = Time of reaction occurred (minutes)

$-r_A$ = Reaction rate (mol/min)

The reactions that occur are:



$-r_A = k \cdot C_A C_B$, since $C_A = C_B$, then

$$-r_A = k C_A^2 = k C^2 \quad (2.11)$$

Eq. (11) $\rightarrow$ Eq. (10)

$$\begin{aligned} \frac{d(V \cdot C)}{dt} &= F_0 \cdot C_0 - V \cdot k \cdot C^2 \\ V \frac{dC}{dt} + C \frac{dV}{dt} &= F_0 \cdot C_0 - V \cdot k \cdot C^2 \end{aligned} \quad (2.12)$$

Eq. (9) $\rightarrow$ Eq. (12)

$$F_0 \cdot t \cdot \frac{dC}{dt} \cdot C = F_0 \cdot C_0 - F_0 \cdot t \cdot k \cdot C^2 \quad (2.13)$$

$$\frac{dC}{dt} = \frac{C_0}{t} - \frac{C}{t} - k \cdot C^2 \quad (2.14)$$

By using the boundary condition at $t=0$, $C=C_0$ and substitution

$U = e^{\int k \cdot C \cdot dt}$ then Eq. (14) becomes :

$$t^2 \frac{d^2 U}{dt^2} + t \frac{dU}{dt} - k \cdot U \cdot C_0 \cdot t = 0 \quad (2.15)$$

Eq. (15) is converted into a Bessel Function by substitution $z = t^{0.5}$, to become:

$$z^2 \frac{d^2 U}{dz^2} + z \frac{dU}{dz} - 4 \cdot k \cdot U \cdot C_0 \cdot z^2 = 0 \quad (2.16)$$

Eq. (16) is a modification of the Bessel Equation which has the following general form:

$$x^2 \frac{d^2 y}{dx^2} + x(a + 2bx^r) \frac{dy}{dx} + [c + dx^{2s} - b(1 - a - r)x \cdot r + b^2 \cdot x^2 \cdot r]. y = 0$$

From the Eq. (2.5) obtained:

$$a = 1$$

$$r = 0$$

$$p = \frac{1}{s} \sqrt{\frac{(1-a)^2}{2}} - c = 0$$

$$b = 0$$

$$s = 0$$

$$p = 0$$

$$c = 0$$

$$d = -4.k.C_0$$

$$\sqrt{\frac{d}{s}} = \text{imaginary}$$

Until the settlement of Eq. (16) are:

$$U = C_1 Z_p (\sqrt{4.k.C_0.z}) + C_2 Z_p (\sqrt{4.k.C_0.z}) \quad (2.18)$$

at $t=0, z=0 \rightarrow Z_p = \sim$

So $C_2 = 0$

$$U = C_1 Z_p (\sqrt{4.k.C_0.z})$$

Because $p = 0$ dan $\sqrt{\frac{d}{s}} = \text{imaginary}$

$$\text{Then } U = C_1 I_0 (\sqrt{4.k.C_0.z})$$

$$\frac{du}{dt} = \frac{d}{dz} C_1 I_0 (\sqrt{4.k.C_0.z}) \quad (2.19)$$

From Sherwood page 178 Eq. (5.83) got

$$\frac{dU}{dt} = C_1 (\sqrt{4.k.C_0.z}) I_0 (\sqrt{4.k.C_0.z}) \quad (2.20)$$

From the original substitution, we get:

$$\frac{dU}{dt} = 2.k.C_2.C_1.I_0 (\sqrt{4.k.C_0.z}) \quad (2.21)$$

Then the Eq. (14) and (15) are obtained:

$$C_1 (\sqrt{4.k.C_0.z}) I_0 (\sqrt{4.k.C_0.z}) = k.C.C_1.I_0 (\sqrt{4.k.C_0.z})$$

$$C = \frac{(\sqrt{4.k.C_0.z}) I_0 (\sqrt{4.k.C_0.z})}{k.C.C_1.I_0 (\sqrt{4.k.C_0.z})}$$

$$C = \frac{C_0 T_1 (2\sqrt{k.C_0 T})}{k.t.T_0 (2\sqrt{k.C_0 T})} \quad (2.22)$$

b. Second Stage

At this stage the process runs continuously, but the steady state has not yet been reached. Can be expressed by:

$$C = f(t) \text{ and } V = \text{constant} \rightarrow dV/dt = 0$$

From the mass balance of the components obtained:

$$\frac{d}{dt} (V.C) = F.C_0 - F.C - k.V.C^2 \quad (2.23)$$

$$V.\frac{dC}{dt} - C.\frac{dV}{dt} = F.C_0 - F.C - k.V.C^2 \quad (2.24)$$

If $T = t - \check{T}$ time, minutes

$$\check{T} = V/F \text{ time constant}$$

Eq. (2.24) will be

$$\frac{dC}{dt} = \frac{C_0}{\check{T}} - \frac{C}{\check{T}} - k \cdot C^2 \quad (2.25)$$

At steady state $C=C_0$

The settlement of the particular Eq (2.25) is $C = C_s$, where C_s is the concentration at steady state.

Substitute $C = C_s + 1/s$

Eq. (2.25) (25) turns into a first-order differential equation which can be solved by the integration factor method.

$$C = C_0 = \frac{1}{B \cdot \exp(AT) - \frac{K}{A}} \quad (2.26)$$

C_1 is the initial concentration of each second stage at the time $t = \check{T}$ obtained by measuring the concentration of the sample.

c. Third phase

At this stage the process is running in a steady state and accumulation = 0 from the component balance, we get:

$$F \cdot C_0 = F \cdot C + V_r \quad (2.27)$$

$$F \cdot C_0 = F \cdot C + V \cdot k \cdot C_s^2 \quad (2.28)$$

$$C_0 = C_s + (V/F) \cdot k \cdot C_s^2 \quad (2.29)$$

$$k \cdot \check{T} \cdot C_s^2 + C_s - C_0 = 0 \quad (2.30)$$


If k is known then C_s can be predicted. On the other hand, if C_s is measured, the value of k can be calculated. Eq. (2.30) is an ordinary algebraic equation and can be solved easily.

2.3 Saponification Reaction

Saponification is generally defined as the hydrolysis reaction of an ester by a strong base, producing a carboxylic acid salt and an alcohol. The term is most commonly associated with soap production, in which triglycerides (fats) react with NaOH or KOH to produce fatty acid salts (soap) and glycerol (Alum, 2024). Nevertheless, the fundamental principle of saponification involves the cleavage of an ester bond by hydroxide ions. The saponification mechanism proceeds through a nucleophilic attack of the hydroxide ion (OH^-) on the carbonyl carbon atom of the ester group, forming a tetrahedral intermediate that subsequently decomposes by releasing an alkoxide group as the leaving group. The strongly basic alkoxide group then abstracts a proton from the carboxylic acid formed, resulting in a carboxylate ion (salt) and an alcohol as the final products (Lestari et al., 2026).

In the Ideal Continuous Flow Reactor laboratory experiment, the saponification reaction investigated is the reaction between ethyl acetate as the

ester source and NaOH as the hydroxide ion source, producing sodium acetate and ethanol. The reaction can be written as follows:



Ethyl acetate and NaOH are selected as reactants because the reaction rate can be readily determined experimentally. The concentration of unreacted NaOH can be directly determined through an acid–base titration using HCl.

In this experiment, the investigation is conducted in two stages: a batch reactor experiment and a continuous reactor (CSTR) experiment. The residual NaOH concentration data obtained from the batch reactor are used to determine the reaction order and reaction rate constant (k) for the saponification of ethyl acetate with NaOH, as described in Sections 2.7 and 2.8. Meanwhile, the residual NaOH concentration data obtained from the continuous reactor experiment are used to validate the mathematical model through the fourth-order Runge–Kutta numerical method, thereby determining the degree of agreement between the experimental results and the model calculations. A thorough understanding of the reaction kinetics enables optimization of operating conditions, such as base concentration, temperature, and reaction residence time, to achieve more efficient reactant conversion (Pasaribu, 2022). In addition, the kinetic data obtained can serve as a basis for designing and controlling the operation of a continuous reactor (CSTR), allowing the process to operate stably and safely while consistently producing products of the desired quality.

2.4 Thermodynamic Overview

Based on reaction : $\text{CH}_3\text{COOC}_2\text{H}_5 + \text{NaOH} \rightarrow \text{CH}_3\text{COONa} + \text{C}_2\text{H}_5\text{OH}$

To determine whether the reaction is exothermic/endothermic, it is necessary to prove it by using the standard heat of formation (ΔH_f°) at 1 atm and 298 K of the reactants and products.

$$\Delta H_{298} = \Delta H_{\text{products}} - \Delta H_{\text{reactants}}$$

The data are known as follows (Smith et al., 2011)

$$\Delta H \text{ CH}_3\text{COOC}_2\text{H}_5 = -445.500 \text{ J/mol}$$

$$\Delta H \text{ NaOH} = -425.609 \text{ J/mol}$$

$$\Delta H \text{ CH}_3\text{COONa} = -726.100 \text{ J/mol}$$

$$\Delta H \text{ C}_2\text{H}_5\text{OH} = -235.100 \text{ J/mol}$$

Therefore :

$$\Delta H_{\text{reaction}} = (\Delta H \text{ CH}_3\text{COONa} + \Delta H \text{ C}_2\text{H}_5\text{OH}) - (\Delta H \text{ CH}_3\text{COOC}_2\text{H}_5 + \Delta H \text{ NaOH})$$

$$= ((-726.100 - 235.100) - (-445.500 - 425.609)) \text{ J/mol}$$

$$= -91.091 \text{ J/mol}$$

Because the ΔH of the reaction is negative, the reaction that takes place is an exothermic reaction that produces heat.



To determine the nature of the reaction whether it goes irreversible or reversible, it can be seen from the value of the reaction equilibrium constant. At room temperature the data obtained (Smith et al., 2011):

$$\Delta G \text{ CH}_3\text{COOC}_2\text{H}_5 = -328.000 \text{ J/mol}$$

$$\Delta G \text{ NaOH} = -379.494 \text{ J/mol}$$

$$\Delta G \text{ CH}_3\text{COONa} = -631.200 \text{ J/mol}$$

$$\Delta G \text{ C}_2\text{H}_5\text{OH} = -168.490 \text{ J/mol}$$

Therefore :

$$\begin{aligned} \Delta G_{\text{reaction}} &= (\Delta G \text{ CH}_3\text{COONa} + \Delta G \text{ C}_2\text{H}_5\text{OH}) - (\Delta G \text{ CH}_3\text{COOC}_2\text{H}_5 + \Delta G \text{ NaOH}) \\ &= (-631.200 - 168.490) - (-328.000 - 379.494) \text{ J/mol} \\ &= -92.196 \text{ J/mol} \end{aligned}$$

$$\frac{d}{dT} \left(\frac{\Delta G}{RT} \right) = \frac{\Delta H}{RT^2}$$

$$\Delta G = RT \ln K$$

$$K \text{ at standard } 298 \text{ K} = e^{(\Delta G/RT)}$$

$$K = e^{\frac{-92.126}{0,314 \cdot 298}} = 1,45 \times 10^{16}$$

From the above data, it can be obtained that the value of the reaction equilibrium constant at a temperature of 298 K is $1,45 \times 10^{16}$. At the operating temperature, the value of K is calculated by the equation:

$$\ln \left(\frac{K}{K'} \right) = \frac{-\Delta H^0}{R} \left(\frac{1}{T} - \frac{1}{T_1} \right)$$

$$T = 27^\circ\text{C (room temperature)} = 300 \text{ K}$$

$$\ln \left(\frac{1,45 \times 10^{16}}{K'} \right) = \frac{-(-91.091)}{8,314} \left(\frac{1}{300} - \frac{1}{298} \right)$$

$$K' = 1,13 \times 10^{16}$$

Because the value of the equilibrium constant is not close to "number" 1, the reaction proceeds in the same direction (irreversible).

2.5 Kinetics Overview

In terms of reaction kinetics, the speed of the saponification reaction of ethyl acetate with NaOH will increase with increasing temperature, stirring and concentration differences. This can be explained by the Arrhenius equation, namely:

$$k = k_0 \cdot e^{-\left(\frac{E_A}{RT}\right)}$$

In this case:

k = Reaction rate constant

k_0 = Pre-exponential factor or frequency

T = Temperature

E_A = Activation energy

R = Ideal gas constant

= 1,98 cal/gm-mol.K

= 1,98 Btu/lb-mol.^oR

= 82,06 cm³. Atm/gm-mol.K

Based on the Arrhenius equation, it can be seen that the reaction rate constant is influenced by the value of the frequency factor or exponential factor, temperature, and activation energy (Levenspiel, 1999).

a. Temperature

Temperature is a factor that affects the reaction rate of the saponification of ethyl acetate with NaOH. The effect of temperature on the saponification reaction rate can be seen in the following figure.




Figure 2.2 Graph of temperature versus reaction rate

(Levenspiel, 1999)

The graph shows the relationship between the reaction rate constant (k) and the reciprocal of absolute temperature ($1/T$) based on the Arrhenius equation, where the reaction rate increases exponentially as temperature increases. Two lines with different slopes represent systems with high and low activation energy (E), indicating that reactions with higher activation energy are more sensitive to changes in temperature.

b. Stirring

Stirring speed is a factor that affects the rate of the saponification reaction of ethyl acetate with NaOH. The effect of stirring speed on the saponification reaction rate can be seen in the following graph.




Figure 2.3 Graph of stirring speed versus reaction rate

(Singh et al., 2016)

The graph illustrates the effect of stirring speed on the reaction rate in catalytic and non-catalytic systems. The reaction rate increases as stirring speed increases up to a certain point, after which it tends to level off. In the catalytic system, the reaction rate is much higher than in the non-catalytic system, indicating the effectiveness of the catalyst in accelerating the reaction. In addition, increasing stirring speed has a more significant effect in the catalytic system, especially at low stirring speeds, suggesting that external diffusion resistance plays an important role.

2.6 Physical and Chemical Properties of Reagents

1. NaOH

Physical Properties

- Molecular weight = 40 g/mol
- Boiling point = 134°C
- Melting point = 318,4°C
- Specific gravity = 2,130 g/mol
- Solubility in 100 parts of cold water 10°C = 42
- Solubility in 100 parts of hot water 100°C = 32

Chemical Properties :

- With Pb(NO₃)₂ it forms a precipitate of Pb(OH)₂ which is soluble in excess reagent, is a strong base, and easily soluble in water.

2. Ethyl Acetate

Physical Properties:

- Specific gravity = 1,356 g/mol
- Boiling point = 85°C
- Molecular weight = 88 g/mol
- Melting point = -111°C

Chemical Properties:

- It is an ester compound formed through the esterification reaction between ethanol and acetic acid, and therefore can undergo hydrolysis reactions under both acidic and basic conditions.
- It is a semi-polar, volatile, non-toxic, and non-hygroscopic solvent.

3. HCl

Physical Properties:

- Atomic mass = 36,45 g/mol
- Density = 3,21 g/ml
- Melting point = -101°C
- Ionization energy = 1250 kJ/mol
- Specific heat = 0,115 kcal/g°C

Chemical Properties:

- HCl produces dense fumes (acid mist) in humid air.
- It is a colorless gas with a sharp and pungent odor.
- It is soluble in alkali hydroxides, chloroform, and ether.
- It is a strong monoprotic acid that completely ionizes in water.
- It has a very high affinity for other elements and can therefore be toxic to the respiratory system.
- It is corrosive and can dissolve various metals, forming chloride salts and hydrogen gas.

2.7 Determining the Order of Reaction

Trial order of reaction in a batch reactor:

1. Given time data (t) and C_A , C_{A0} is C_A at $t=0$
2. Create data $-\ln (C_A/C_{A0})$
3. First guess the 'first order reaction' by graphing $-\ln (C_A/C_{A0})$ vs t, the resulting graph should be straight.

- If the result of the graph is not straight then guess the 'second order reaction' from the graph between $1/C_A$ vs t , the result of the graph must be straight (If $C_{A0} = C_{B0}$)
 - If the graph results are not straight then guessing the 'second order reaction' from the graph between $\ln C_B/C_A$ vs t , the graph results must be straight.
- Form the equation $y = a + bx$, where a = intercept and b = slope of the $\log t$ vs $\ln C_{A0}$ graph.




Figure 2.2 1st Order reaction trial graph




Figure 2.3 2nd Order reaction trial graph ($C_{A0} = C_{B0}$)




Figure 2.4 2nd Order reaction trial graph ($C_{A0} \neq C_{B0}$)





Figure 2.5 n Order reaction trial graph

Process Laboratory

2.8 Calculating the Saponification Reaction Constant (k) Ethyl Acetate with NaOH



1st Order Reaction

$$-r_a = -\frac{dC_a}{dt} = k \cdot C_a$$

$$\int_{C_{a0}}^{C_a} -\frac{dC_a}{C_a} = \int_0^t k \cdot dt$$

$$-\ln[C_a]_{C_{a0}}^{C_a} = k \cdot t$$

$$-(\ln C_a - \ln C_{a0}) = k \cdot t$$

$$-\ln \frac{C_a}{C_{a0}} = k \cdot t$$

$$y = mx$$

2nd Order Reaction

The reaction rate equation if $C_a = C_b$:

$$-r_a = -\frac{dC_a}{dt} = k \cdot C_a \cdot C_b$$

$$-\frac{dC_a}{dt} = k \cdot C_a^2$$

$$-\frac{dC_a}{dt} = k \cdot dt$$

$$\int_{C_{a0}}^{C_a} -\frac{dC_a}{C_a^2} = \int_0^t k \cdot dt$$

$$\left[\frac{1}{C_a} \right]_{C_{a0}}^{C_a} = k \cdot t$$

$$\frac{1}{C_a} - \frac{1}{C_{a0}} = k \cdot t$$

$$\frac{1}{C_a} = k \cdot t + \frac{1}{C_{a0}}$$

$$y = mx + c$$

The reaction rate equation if $C_a \neq C_b$:

$$r_a = -\frac{dC_a}{dt} = k \cdot C_a \cdot C_b$$

$$-r_a = -C_{a0} \frac{dX_a}{dt} = k(C_{a0} - C_{a0} \cdot X_a)(C_{b0} - C_{b0} \cdot X_a)$$

$$\frac{dX_a}{(1-X_a)(M-X_a)} = C_{a0} \cdot k \cdot t, \text{ in this case } M = \frac{C_{b0}}{C_{a0}}$$

$$\frac{1}{C_{a0}(M-1)} \int_0^{X_a} \frac{dX_a}{1-X_a} - \int_0^{X_a} \frac{dX_a}{M-X_a} = k \cdot t$$

$$\ln \frac{M-X_a}{M(1-X_a)} = C_{a0} \cdot (M-1)k \cdot t \text{ or } \ln \frac{M-X_a}{M(1-X_a)} = (C_{a0} - C_{b0})k \cdot t$$

$$\ln \frac{M-X_a}{M(1-X_a)} = \ln \frac{C_b - C_{a0}X_a}{C_{b0}(1-X_a)}$$

$$\ln \frac{M-X_a}{M(1-X_a)} = \ln \frac{C_b C_{a0}}{C_{b0}(1-X_a)C_{a0}}$$

$$\ln \frac{M-X_a}{M(1-X_a)} = \ln \frac{C_b C_{a0}}{C_{b0} C_{a0}} = \ln \frac{C_b}{M \cdot C_a}$$


$$\ln \frac{C_b}{M \cdot C_a} = (C_{b0} - C_{a0})k \cdot t$$

$$\ln \frac{C_b}{C_a} = (C_{bo} - C_{ao})k \cdot t + \ln M$$

$$y = mx + c$$

The value of k is obtained from the least square method. Where the value of k is the value of m.

(Levenspiel, 1999)



CHAPTER III

PRACTICUM METHOD

3.1 Materials and Tools Used

3.1.1 Materials and Tools Used

1. NaOH
2. Ethyl acetate
3. HCl
4. MO Indicator
5. Aquadest

3.1.2 Tools used

1. Pipette
2. Batch reactor
3. Measuring cup
4. Burette
5. Statives and Clamps
6. Erlenmeyer
7. Continuous flow reactor circuit

3.2 Experiment Circuit Drawing

a. Batch Process




Figure 3.1 Batch process main tool drawing

Information:

1. batch reactor
2. *Stirrer*
3. Stative

b. Continuous Process




Figure 3.2 Continuous process main tool drawing

Information :

1. Continuous reactor
2. Stirrer
3. Statives
4. Reactor tank

3.4 Response Test Results

The residual NaOH concentration that can be observed with the titrant concentration of HCl to ET (Endpoint Titration)

3.5 Trial Procedure


a. Batch trial

1. Prepare the required reagents: ethyl acetate at ... N, HCl at ... N, and NaOH at ... N.
2. Add 1 L each of ... N ethyl acetate and ... N NaOH.
3. Take a 5 mL sample every ... minutes, then add 3 drops of methyl orange (MO) indicator to the sample and titrate with ... N HCl until a reddish-orange color is obtained. The titration is stopped when the volume of titrant used remains constant for three consecutive measurements.
4. The value of C_A (residual NaOH concentration) can be determined through calculation.
5. Repeat Steps 1–4 using different variables.

a. Continuous Trial

1. Prepare the required reagents: ethyl acetate at ... N, HCl at ... N, and NaOH at ... N.
2. Add ethyl acetate and NaOH into their respective feed tanks.

3. Adjust the valves of each reactant to feed them into the empty CSTR and maintain a constant flow rate for both reactants.
4. Take a 5 mL sample every ... minutes, then add 3 drops of methyl orange (MO) indicator to the sample and titrate with ... N HCl until a reddish-orange color is obtained. The titration is stopped when the volume of titrant used remains constant for three consecutive measurements.
5. The value of CA (residual NaOH concentration) can be determined through calculation.
6. Repeat Steps 1–5 under moderate and high stirring conditions.



Process

Laboratory

BIBLIOGRAPHY

- Alum, B. N. (2024). Saponification process and soap chemistry. *INOSR Applied Sciences*, 12(2), 51-56.
- Lestari, D. I., Sari, D. A., Putri, Z. P., & Maranggi, I. U. (2026). Kinetics of irreversible saponification reactions in terms of changes in reactant concentration: kinetika reaksi saponifikasi irreversible ditinjau dari perubahan konsentrasi reaktan. *Jurnal Kimia dan Rekayasa*, 6(2), 44-51.
- Levenspiel, O., 1999. *Chemical reaction Engineering 3rd ed*, Mc. Graw Hill Book Kogakusha Ltd, Tokyo.
- Pasaribu, M. H. (2022). Studi kinetika kimia pada reaksi saponifikasi etil asetat dengan metode titrasi dan konduktometri: chemical kinetic study on saponification reaction of ethyl acetate using titration and conductometric methods. *Jurnal Jejaring Matematika dan Sains*, 4(2), 58-64.
- Singh, G., Nakade, P. G., Mishra, P., Jha, P., Sen, S., & Mondal, U. (2016). Kinetic investigation on liquid-liquid-solid phase transfer catalyzed synthesis of dibenzyl disulfide with H₂S-laden monoethanolamine. *Journal of Molecular Catalysis A: Chemical*, 78 - 86.
- Smith, J. M. ,VanNess, H. C., Abbott, M. M. 2011. *Introduction to Chemical Engineering Thermodynamics*, 6th ed., McGraw-HillCo., Singapore

Process

Laboratory

HAZARD IDENTIFICATION AND RISKS ANALYST
MATERIAL: THE IDEAL CONTINUOUS FLOW REACTOR

HAZARD IDENTIFICATION (HI)											
A	Mechanic		D	Environment		E	Chemical Material		G	Other hazard	
A1	Manual handling	√	D1	Noise		E1	Toxicity	√	G1	Compressed gas	
A2	Moving part		D2	Vibration	√	E2	Irritant		G2	Ionic radiation	
A3	Rotating part	√	D3	Light		E3	Corosive	√	G3	UV radiation	
A4	Cutting		D4	Humidity		E4	Carcinogenic		G4	Fatigue	
B	Biology		D5	Temperatuer	√	E5	Flammable	√	G5	Narrow space	
B1	Bacterial		D6	Travel hazard		E6	Explosive		G6	Crowded	
B2	Virus		D7	Slippery surface		E7	Cryogenics		G7	Thermometer	
B3	Fungi		D8	Solid waste		F	Equipments				
C	Electricity		D9	Air qualiy		F1	Pressure vessel				
C1	High voltage		D10	Solitary work		F2	Hot equipments				
C2	Electrostatic		D11	Splash/water drop/flood	√	F3	Laser				
C3	Cable	√	D12	Powder spill		F4	Glass vessel				

RISKS							
IB	Risks (after action)				Risk Identification	Action to Minimize Risks	First Aid Action
	High	Medium	Low	Minimal			
1. PREPARATION							
			√		NaOH, HCl, and Ethyl Acetate preparation	- Make sure to always wear gloves when taking solutions, especially HCl solutions in the acid room - Always turn on the fan in the acid chamber, if you are taking HCl solution - Always wear complete lab equipment, such as lab coats, gloves, shoes, and medical masks and a mask	- Wash hands if exposed to related solutions - If it accidentally hurts, immediately ask for further help
2. MAIN EXPERIMENTS							
			√		When starting the pump there is a risk of electric shock, or the pump will catch fire	- Make sure the cable / plug is not wet - Always wear complete lab equipment, such as lab coats, gloves, shoes, and medical masks	- Turn off the power source - Push the victim's body with an insulating object - Immediately ask for first aid
3. ANALYSIS/FINAL STAGE							

Process

Laboratory

					NaOH, HCl, and Ethyl Acetate waste disposal	- Make sure the cable / plug is not wet - - Always make sure each of chemicals have their own place to be disposed - Always wear complete lab equipment, such as lab coats, gloves, shoes, and medical masks - Always wear complete lab equipment, such as lab coats, gloves, shoes, and medical masks	- Wash hands if exposed to related solutions - If it accidental hurts, immediately ask for further help
--	--	--	--	--	--	---	--

Process

Laboratory